
Internet: www.basko.com
 E-mail: sales@basko.com
 Fax: +31 (0) 75 - 612 63 73
 Phone: +31 (0) 75 - 613 15 13
 NL-1505 HX Zaandam
 Pieter Liefincxweg 16

Basko Healthcare


basko.com

Basko Healthcare

FreeMotion

Alignment Unit

AK-104

- NL** Gebruiksaanwijzing
- DE** Gebrauchsanleitung
- ENG** Instructions for Use
- FR** Utilisation et application
- ES** Uso y aplicación
- IT** Usi e applicazioni
- SK** Bruksanvisning
- DK** Bruksanvisning
- NO** Bruksanvisning
- FI** Sovitusohjeet
- CZ** Návod k použití
- SI** Navodila za uporabo

Basko Healthcare
basko.com

NL

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ENG

Application Guidelines
 Proper construction with associated joint alignment is critical to a properly functioning orthosis. This alignment kit facilitates joint alignment during the production of an orthosis with the FreeMotion joint. If applicable, also use in combination with a MultiMotion hinge for which we recommend the additional AK-102 alignment kit or MultiStatic hinge with alignment kit AK-103.

The certified technicians producing the orthosis are responsible for the production technique chosen and the correct functioning of the orthosis. Our advice is to use the product-specific alignment unit at all times during production.

Prior to production, all possible openings on the joint should be sealed with suitable material to prevent the joint from filling up with resin or similar.

Additional product information
 The alignment unit is suitable for multiple use.

Maintenance instructions
 After use, inspect that all parts are undamaged, clean and free of plaster and similar.

Instructions
 The alignment unit for the FreeMotion joint consists of 2 separate parts:

# AK-100 Basic T-Bar	# AK-104 FM kit
First check that all parts shown in Figure 1 are complete and undamaged in the package:	

Package contents AK-100:
A) 1 x T-bar

Package contents AK-104:
A) 2 x long bar profile
B) 2 x short bar profile
C) 2 x spacer
D) 4 x cylindrical screw with hexagon socket
E) 1 x Allen key (3 mm)
F) Instruction for use

Step 1: While casting the patient, mark the joint axes. After removing and resealing the negative, perforate the axis markings at the desired position and then insert the T-bar through the perforation. After piercing, check that the axis is horizontal and close the negative around the T-bar (Fig. 2).

Step 2: Fill the negative with plaster possibly in combination with Vermiculite (art. no. 17) and let it harden well. After curing, remove the plaster negative and the T-bar. Then correct the plaster negative as needed and allow it to dry (Fig. 3).

Step 3: Insert the cylinder head screw through the opening of the FreeMotion joint axis (Fig. 4).

Step 4: Then (depending on the M-L width of the plaster) screw the short or long bar profile onto the cylinder head screw with the supplied Allen wrench. Use the spacer as needed depending on the contour of the plaster model and desired joint offset (Fig. 5).

Step 5: Using the bar profile, place the joint parallel to the "joint" in the square T-bar opening (Fig. 6).

Step 6: Where necessary, bend the bars to the contours of the plaster positive (Fig. 7).

For bilateral construction of the Orthosis, repeat steps 3 through 6 on the opposite side. After assembling the Orthosis, unscrew the cylinder head screw and remove the bar profile.

The alignment unit is suitable for multiple use.

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