# BINKY

Universal Free Motion Joint

Professional instructions

**Basko**  
Healthcare

[basko.com](http://basko.com)

|                                     |    |
|-------------------------------------|----|
| Technische gebruiksaanwijzing ..... | 4  |
| Technische Gebrauchsanleitung ..... | 6  |
| Technical manual .....              | 8  |
| Manuel technique .....              | 10 |
| Manual técnico .....                | 12 |
| Manuale tecnico .....               | 14 |
| Teknisk manual .....                | 16 |
| Teknisk manual .....                | 18 |
| Teknisk manual .....                | 20 |
| Tekninen käsikirja .....            | 22 |

NL

DE

ENG

FR

ES

IT

SE

DK

NO

FI

## BINKY

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## BINKY

Please read these instructions carefully before using the product.  
Pay particular attention that the safety instructions are followed.

### Usage guidelines

- The product may only be fitted by or under the supervision of a competent professional
- Special precautions must always be taken for patients with a visual impairment, cognitive limitations and/or reduced sensibility in the upper and/or lower extremity
- Changes, modifications and adjustments to the product that are not described in these instructions, are done at the sole responsibility of the person carrying them out

### Intended use / Application

These universal joints without flexion and / or extension stops are used pairwise for the production of orthoses made out of thermoplastics or separately as a free motion joint. BINKY joints are suitable for use by both children and adults in orthoses for wrist, elbow, and ankle joint. The physician is the one who determines the indications and also determine the duration of treatment.

### Indications

- Paralytic syndromes (other) (ICD 10: G83.-).

### Contraindications:

There are no known contraindications.

### General safety instructions

*Risk of injury due to excessive load on the joints:* structural limits and specific values i.e. body weight may determine size selection of the joints. Excessive flexion loads may lead to breakage of the joints.

*Risk of injury due to joint alignment deviations:* during fabrication, correct joint alignment must be assured in the frontal and sagittal plane to ensure smooth movement of the joints.

*Risk of injury by using the wrong materials:* the joints have been developed for usage in orthoses made of thermoplastic material. The use of the joint in orthotics made of other materials is not recommended.

*Risk of injury as a result of incorrect use:* the professional should verify correct functioning of the joints prior to delivery of the orthosis to the user. The patient and his/her care giver(s) should also be instructed on donning and doffing. Any settings may only be changed by the technician, therapist or doctor.

### Used materials

- Stainless steel

### Joint selection

BINKY joints are available in 2 sizes. Directive is to apply the 814101 when used on small joints and the 814102 for larger joints. Other than that, it is up to the professional judgment of the technician to select the desirable joint based on knowledge and experience.

## Size chart

| Art. Nr. | Unit   | Length | Width | Thickness | Body weight |
|----------|--------|--------|-------|-----------|-------------|
| 814101   | 1 pair | 62 mm  | 12 mm | 2,5 mm    | < 35 kg     |
| 814102   | 1 pair | 76 mm  | 18 mm | 2,5 mm    | < 100 kg    |

## Maintenance recommendation

We recommend a service check on functionality, wear and tear once every six months.

## Cleaning instructions

The joint may be cleaned with naphtha or similar. Try to keep the joint free from dust and dirt as much as possible.

## Package content

- 1 set BINKY joints
- Technical instructions

Check the package label to make sure the article number corresponds with the content of the package.

## Handling

- All general applicable technical and anatomical directives should be taken into account when manufacturing the plaster model, determination of pivot points, alignment, assembly and orthosis construction
- Each joint needs to be mounted onto the orthosis at all 4 points

## Liability

Basko Healthcare as the manufacturer is solely liable if the product is used in accordance with the prescribed conditions and the purpose for which it is intended. Basko Healthcare recommends handling the product in accordance with applicable rules and performing maintenance as described in these user instructions. Non-compliance with the above may cause the BINKY joints to not function or function inadequately and may impact guarantee.

## CE-conformity

These products meet requirements in Directive 93/42/EEC for medical devices. These products have been included in Class I based on classification criteria for medical devices as per annex IX of the Directive. Basko Healthcare therefore wrote up this conformity statement entirely at its own responsibility as per annex VII of the Directive.

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