
Internet: www.basko.com
 E-mail: sales@basko.com
 Fax: +31 (0) 75 - 612 63 73
 Phone: +31 (0) 75 - 613 15 13
 NL-1505 HX Zaandam
 Pieter Liefinckweg 16

Basko HealthCare


basko.com

Basko HealthCare

Swing Phase Lock 2

Alignment Unit

AK-101

- NL Gebruiksaanwijzing
- DE Gebrauchsanleitung
- ENG Instructions for Use
- FR Mode d'emploi
- ES Uso y aplicación
- IT Usi e applicazioni
- SK Bruksanvisning
- DK Bruksanvisning
- NO Bruksanvisning
- FI Sovitusohjeet
- CZ Návod k použití
- SI Navodila za uporabo

Basko HealthCare
 basko.com

NL

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.....

ENG

Application Guidelines
 Proper construction with associated joint alignment is critical to a properly functioning orthosis. This alignment kit facilitates joint alignment during the production of an orthosis with Swing Phase Lock 2 (SPL 2) joint. This alignment kit can also be used in conjunction with a Swing Phase Control (SPC) joint, if applicable.

The certified technicians producing the orthosis are responsible for the production technique chosen and the correct functioning of the orthosis. Our advice is to use the product-specific alignment unit at all times during production.

Prior to production, all possible openings on the joint should be sealed with suitable material to prevent the joint from filling up with resin or similar.

Additional product information
 The alignment unit is suitable for multiple use. Two types of joint dummies are separately available: a steel dummy set 029860 and plastic dummies IQ200R-D or IQ200L-D.

Maintenance instructions
 After use, inspect that all parts are undamaged, clean and free of plaster and similar.

Instructions
 The alignment unit for the Swing Phase Lock 2 hinge consists of 2 separate parts:

AK-100 Basic T-Bar # AK-101 SPL2 kit

First check that all parts shown in Figure 1 are complete and undamaged in the package:

Package contents AK-100:

A) 1 x T-bar

Package contents AK-101:

- A) 2 x threaded bar profile
- B) 4 x cylindrical screw with hexagon socket
- C) 1 x Function unit dummy
- D) 1 x Allen key (3mm)
- E) 1 x Allen key (6 mm)
- F) Instruction for use

Step 1: While casting the patient, mark the joint axes. After removing and resealing the negative, perforate the axis markings at the desired position and then insert the T-bar through the perforation. After piercing, check that the axis is horizontal and close the negative around the T-bar (Fig. 2).

Step 2: Fill the negative with plaster possibly in combination with Vermiculite (art. no. 17) and let it harden well. After curing, remove the plaster negative and the T-bar. Then correct the plaster negative as needed and allow it to dry (Fig. 3).

Step 3: Use the Allen wrench supplied with the SPL2 joint to remove the blue cover plate on the front and the stainless steel cover plate on the back of the joint (Fig. 4).

Step 4: Next, loosen the center axle screw from the axis using the supplied 3mm Allen wrench and push the axis out of the joint housing (Fig. 5). Then remove the satellite and function unit from out of the joint as described in the SPL2 joint instructions for use. After removing, insert the function unit dummy.

Step 5: Next, insert a hexagonal cylindrical screw through the central axis on the front of the joint and screw it in place with the supplied 6 mm Allen wrench provided onto the bar profile at the rear of the joint. (Fig. 6).

Step 6: Next, install the uprights (Art. No. IQ023516 / IQ023616) or stirrups (Art. No. IQ161 / IQ162) distal and proximal onto the SPL2 joint (Fig. 7).

Step 7: Then laterally place the SPL2 joint parallel to the "joint" using the bar profile in the opening where the T-bar has been (Fig. 8).

Step 8: Where necessary, bend the uprights/stirrups to the follow the contours of the plaster (Fig. 9).

For a bilateral construction (with SPC joint) of the Orthosis, repeat steps 4 through 8 on the medial side noting that the SPC joint has only a front cover.

After the orthosis is fabricated, remove the alignment unit following steps 4 through 6 in reverse order and follow the instructions in the SPL2 Knee Joint User Manual to adjust and check proper function.

The alignment unit is suitable for multiple use.

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