

Basko
Healthcare



Clubycare Boots



basko.com

Clubycare Boots

Clubycare boots are made of soft Nappa leather and offer high comfort, tolerance and user compliance. These boots are specially designed for the treatment method according to Dr. Ponseti. The boots are mounted on to the foot plates of the P.E.A.care splints.

Features

- Available with lace- or velcro closure
- Padded and increased calcaneal cupping
- Open heel construction
- Open toe construction
- Sold **per piece**
- Elongated heel strap (stitched markers, unstretchable)
- Works perfectly in conjunction with P.E.A.care splints
- Material: nappa leather; sole: rubber; lace: cotton; padding: EVA-foam

Indications

- Arthrogryposis multiplex congenita (ICD10: Q74.3)
- Pes equinovarus congenitus (ICD10: Q66.0)



Sizing chart - Clubycare Boots

Art. No. without abduction left with laces	Art. No. without abduction right with laces	Art. No. without abduction left with velcro	Art. No. without abduction right with velcro	Art. no. with abduction left with Velcro	Art. no. with abduction right with Velcro	Shoe size	Shoe Length	Foot width
BA905L-12	BA905R-12	BA908L-12	BA908R-12	-	-	12	72 mm	40 mm
BA905L-13	BA905R-13	BA908L-13	BA908R-13	-	-	13	79 mm	42 mm
BA905L-14	BA905R-14	BA908L-14	BA908R-14	BA705L-14	BA705R-14	14	85 mm	43 mm
BA905L-15	BA905R-15	BA908L-15	BA908R-15	BA705L-15	BA705R-15	15	92 mm	45 mm
BA905L-16	BA905R-16	BA908L-16	BA908R-16	BA705L-16	BA705R-16	16	98 mm	46 mm
BA905L-17	BA905R-17	BA908L-17	BA908R-17	BA705L-17	BA705R-17	17	107 mm	47 mm
BA905L-18	BA905R-18	BA908L-18	BA908R-18	BA705L-18	BA705R-18	18	113 mm	48 mm
BA905L-19	BA905R-19	BA908L-19	BA908R-19	BA705L-19	BA705R-19	19	120 mm	50 mm
BA905L-20	BA905R-20	BA908L-20	BA908R-20	BA705L-20	BA705R-20	20	126 mm	51 mm
BA905L-21	BA905R-21	BA908L-21	BA908R-21	BA705L-21	BA705R-21	21	133 mm	52 mm
BA905L-22	BA905R-22	BA908L-22	BA908R-22	BA705L-22	BA705R-22	22	139 mm	54 mm
BA905L-23	BA905R-23	BA908L-23	BA908R-23	BA705L-23	BA705R-23	23	146 mm	55 mm
BA905L-24	BA905R-24	BA908L-24	BA908R-24	BA705L-24	BA705R-24	24	152 mm	56 mm
BA905L-25	BA905R-25	BA908L-25	BA908R-25	BA705L-25	BA705R-25	25	159 mm	58 mm
BA905L-26	BA905R-26	BA908L-26	BA908R-26	BA705L-26	BA705R-26	26	166 mm	59 mm
BA905L-27	BA905R-27	BA908L-27	BA908R-27	BA705L-27	BA705R-27	27	172 mm	60 mm
BA905L-28	BA905R-28	BA908L-28	BA908R-28	BA705L-28	BA705R-28	28	179 mm	62 mm
BA905L-29	BA905R-29	BA908L-29	BA908R-29	BA705L-29	BA705R-29	29	185 mm	63 mm
BA905L-30	BA905R-30	BA908L-30	BA908R-30	-	-	30	192 mm	64 mm

Sold per piece. An essential function of the Clubycare shoe is the heel strap, which ensures the heel is correctly positioned against the heel counter of the shoe. Therefore, sufficient growth reserve should be taken into account (max. 15 mm).