



Rhizocare - Thumb Splint

Rhizocare - Thumb Splint

This thumb splint is used both conservative as well as post-operative and intended to immobilize the wrist, thumb, saddle and thumb basal joint.

Indications

- Arthritis, rheumatoid
- Distortion
- Ligament Rupture (Skier's Thumb)
- Radialis paresis
- Soft Tissue irritation

Features

- Anatomically shaped
- Integrated Spring Steel reinforcement
- Adjustable Thumb strap
- Material: Polyester, Nylon, Spring Steel



Sizing Chart

Item number		Size	Length	Wrist circumference	
Right	Left				
8711R-S	8711L-S	small	16,5 cm	12 - 15 cm	
8711R-M	8711L-M	medium	16,5 cm	15 - 17 cm	
8711R-L	8711L-L	large	16,5 cm	17 - 19 cm	
8711R-XL	8711L-XL	X-large	16,5 cm	19 - 21 cm	

Measurement: circular at the widest part of the wrist