

Mendivil Boots

Mendivil boots are made of soft leather and offer high comfort, tolerance and user compliance. The soles of the boots are easy to screw on to the foot plates of the P.E.A.care splint. These screws are supplied with the splint. Mendivil boots come in two types - with lace closure* and with velcro closure** - in different sizes.

* with or without abduction.
** without abduction.



Features

- Reversible padded tongue
- Padded and increased calcaneal cupping
- Open heel construction
- Open toe construction
- Elongated fastening strap with clearly stiched markers
- Works perfectly in conjunction with P.E.A.care splints
- Sold in pairs
- Material: leather

Size chart - Mendivil Boots

Article number lace closure with abduction	Article number lace closure without abduction	Article number velcro closure without abduction	Shoe size	Shoe length	Joint width M. 1-5 with abduction (ME705)	Joint width M. 1-5 without abduction (ME905/ME908)
-	ME905-12	ME908-12	12	72 mm	-	40 mm
ME705-13	ME905-13	ME908-13	13	79 mm	46 mm	42 mm
ME705-14	ME905-14	ME908-14	14	85 mm	47 mm	43 mm
ME705-15	ME905-15	ME908-15	15	92 mm	49 mm	45 mm
ME705-16	ME905-16	ME908-16	16	98 mm	50 mm	46 mm
ME705-17	ME905-17	ME908-17	17	107 mm	51 mm	47 mm
ME705-18	ME905-18	ME908-18	18	113 mm	52 mm	48 mm
ME705-19	ME905-19	ME908-19	19	120 mm	54 mm	50 mm
ME705-20	ME905-20	ME908-20	20	126 mm	55 mm	51 mm
ME705-21	ME905-21	ME908-21	21	133 mm	56 mm	52 mm
ME705-22	ME905-22	ME908-22	22	139 mm	58 mm	54 mm
ME705-23	ME905-23	ME908-23	23	146 mm	59 mm	55 mm
ME705-24	ME905-24	ME908-24	24	152 mm	60 mm	56 mm
ME705-25	ME905-25	ME908-25	25	159 mm	62 mm	58 mm
ME705-26	ME905-26	ME908-26	26	166 mm	63 mm	59 mm
ME705-27	ME905-27	ME908-27	27	172 mm	64 mm	60 mm
ME705-28	ME905-28	ME908-28	28	179 mm	66 mm	62 mm
ME705-29	ME905-29	ME908-29	29	185 mm	67 mm	63 mm
-	ME905-30	ME908-30	30	192 mm	68 mm	64 mm
-	ME905-31	ME908-31	31	198 mm	70 mm	66 mm
-	ME905-32	ME908-32	32	205 mm	71 mm	67 mm
-	ME905-33	ME908-33	33	212 mm	72 mm	68 mm

Basko
Healthcare

Pediatric Clubfoot Treatment



P.E.A.care Splint and
Mendivil Boots

Basko Healthcare
Pieter Liefinckweg 16 | 1505 HX Zaandam
Tel.: +31 (0) 75 - 613 15 13 | Fax: +31 (0) 75 - 612 63 73
E-mail: verkoop@basko.com | Internet: www.basko.com

basko.com

8500703 / Version: 12/2013

basko.com

For the treatment of clubfoot with babies and toddlers Basko Healthcare provides orthotic components such as the P.E.A.care splint and Mendivil boots. This combination in particular is used in the treatment of clubfoot according to the Dr. Ponseti method.

Indication

- Pes equinovarus (Clubfoot) (ICD10: Q66.8/ M21.57)

P.E.A.care Splint

The special features, it's design, and particularly the lightweight of the P.E.A.care splints are just a few of the versatile arguments for the treatment with one of these components. Three different models are available in multiple sizes:

Model "Static"

Features

- Adjustable splint width
- Adjustable rotation angle setting (15° to 85°)
- Dorsalflexion with "P.E.A.care-extension setting" possible
- Optional Quick Release conversion-kit available
- Combinable with Mendivil boots
- Weight: 113 - 124 gram
- Material: aluminum



Model "Static Quick Release"

Features

- Adjustable splint width
- Adjustable rotation angle setting (15° to 85°)
- Dorsalflexion with "P.E.A.care-extension setting" possible
- Quick Release mechanism provides easy donning and doffing of the boots
- Combinable with Mendivil boots
- Weight: 210 - 220 gram
- Material: aluminum



Model "Dynamic"

Features

- Adjustable splint width
- Adjustable rotation angle setting (15° to 85°)
- Dorsalflexion with "P.E.A.care-extension setting" possible
- Optional Quick Release conversion-kit available
- Combinable with Mendivil boots
- Weight: 240 - 250 gram
- Material: aluminum



Model "Dynamic Quick Release"

Features

- Adjustable splint width
- Adjustable rotation angle setting (15° to 85°)
- Quick Release mechanism provides easy donning and doffing of the boots
- Dynamic joint for more freedom of movement (treading or crawling)
- Combinable with Mendivil boots
- Weight: 250 - 270 gram
- Material: aluminum



Size chart

Art. no. "Static"	Art. no. "Static Quick Release"	Art. no. "Dynamic"	Art. no. "Dynamic Quick Release"	Size	Splint width	Shoulder width
2800-1	2810-1	2830-1	2840-1	1	220 - 290 mm	185 - 260 mm
2800-2	2810-2	2830-2	2840-2	2	270 - 400 mm	245 - 370 mm

Measurement: The appropriate size is determined by the shoulder width.

Orderinfo accessories

Quick Release-kit (Only Static)	Extension setting
Art. no. 2819*	Art. no. 2828*

*Set of 2 pieces.